

Exhibit A

11/15/95

Minipreps of DIRECTED DNAs in CONTIG.

1) 205-A2, 2) 205-A4, 3) 205A7 4) BTS1-A6. 5) BTS3A-A3.
6) BTS3A-D6 7) BTS3A-E8 8) BTS3B-E12 9) BTS3B-F4.

- do 3.0 ml minipreps of each. \Rightarrow 30 μ l H₂O R/Vase
2 μ l DNA
1.5 μ l 10x H
1.0 μ l EcoRI (10u)
0.5 μ l H₂O.

