

Exhibit A

11/15/95

Minipreps of DIRECTED DNAs in CONTIG.

- 1) 205-A2, 2) 205-A4, 3) 205A7 4) BT51-A6. 5) BT53A-A3.
6) BT53A-D6 7) BT53A-E8 8) BT53B-E12 9) BT53B-F4.

- do 3.0 ml minipreps of each. \Rightarrow 30 μ l H₂O R/Vase
2 μ l DNA
1.5 μ l 10X H
1.0 μ l EcoRI (10u)
0.5 μ l H₂O.

